

**Sealapex Xpress Recall (8/10)**

Sybron Dental Specialties recently initiated a voluntary recall on Sealapex Xpress (Part Number 33639 and Lot Numbers: 6-1314, 7-1221, and 7-1258). Sealapex Xpress is intended as a root canal filling material that is used in conjunction with gutta percha or silver endodontic points. The company is taking this action due to the presence of crystals in the catalyst. This crystallization may affect proper canal filling by preventing gutta percha points from reaching the bottom of the root canal which may lead to an inadequate seal. Clinics should return the affected products to SybronEndo; any product returned will be replaced at no charge. For additional questions about this recall, please contact SybronEndo Customer Care Department at (800) 346-3838.